


FITTING PROGRESSIVE LENSES IN A  
SUNWEAR WRAP FRAME IS  
**NO LONGER  
A PROBLEM**

DECENTERED POLY IMAGE® NUPOLAR® GRAY OFFERS  
LARGER EFFECTIVE LENS SIZE IN CRITICAL AREAS



**IMAGE® WRAP  
DECENTERED DESIGN**

STANDARD PROGRESSIVE DESIGN DOES NOT OFFER ENOUGH LENS  
CUT OUT MATERIAL TO FIT MANY WRAPPED SUNWEAR FRAMES




IMAGE®

NUPOLAR®  
*polarized lenses*

## Benefits of Decentered design

- Available in +8 BC for "Wrap Frames"
- Decentered 7mm to allow for maximum cut-out opportunity
- An effective blank size of 83mm
- No secondary calculations necessary for use
- Great for narrow PD's
- Now the ability to match premium Image® NuPolar® Poly technology with today's premium sunglass frames
- No need to compromise technology and performance on Rx sunglasses due to cut-out issues



## Availability

IMAGE®

NUPOLAR®  
polarized lenses

IMAGE® PROGRESSIVE POLYCARBONATE - GRAY 3 - HARD COATED									
diameter	nominal curve	true curve	back curve	edge	center	MRP dec.	seg inset	fitting cross above MRP	add range
76	8.5	8.31	9.5	10.5	9	7	2.5	2	1.00-3.00